


Nitinol Stone Extraction Baskets

Single Use

Niti-8-Extractor

- Kink Resistant Nitinol Baskets
Reliable and controlled radial deployment in the bile duct
- Rotatable Basket in 8-Wire Design Facilitated by Synchronous Basket / Handle Rotation Effect
Secure capture and optimal retrieval of small stones and fragments
- Double Lumen Guide Wire Compatible Versions
Providing excellent ductal access
Precise maneuverability through & positioning in challenging bile duct anatomy



Nitinol Stone Extraction Baskets

Single Use

Features & Benefits

Kink Resistant Nitinol Baskets

Reliable and controlled radial deployment in the bile duct

Rotatable Basket in 8-Wire Design Facilitated by Synchronous Basket / Handle Rotation Effect
Secure capture and optimal retrieval of small stones and fragments

Double Lumen Guide Wire Compatible Versions
Providing excellent ductal access

Precise maneuverability through & positioning in challenging bile duct anatomy

Ergonomic Handle Design
Safe basket manipulation

45° Injection Port & Integrated Sealing
User-friendly injection of contrast medium

Double Lumen Guide Wire Compatible Stone Extraction Basket

Niti-8-Extractor – Nitinol Stone Extraction Basket, Working Length 180 cm

Catalog No.	Description	Basket Diameter	Basket Length	Outer Diameter	Max. Accepted Guide Wire
GEB-14-27-820	8-Wire Basket	20 mm	34 mm	2.7 mm	.035"
GEB-14-27-830	8-Wire Basket	30 mm	49 mm	2.7 mm	.035"

Single Lumen Stone Extraction Basket

Niti-8-Extractor – Nitinol Stone Extraction Basket, Working Length 180 cm

Catalog No.	Description	Basket Diameter	Basket Length	Outer Diameter	Max. Accepted Guide Wire
GEB-13-19-820	8-Wire Basket	20 mm	34 mm	1.9 mm	-
GEB-13-19-830	8-Wire Basket	30 mm	49 mm	1.9 mm	-

Caution: This device is not detachable and (must not be used) compatible with the StoneBuster® Lithotripter or any other mechanical lithotripsy system.

Distributed by:

Comercializadora Endosistemas Ltda.
Suecia N° 3243, Ñuñoa, Stgo.
Fonos: 222090262 - 227855826
ventas@endosistemas.cl

Manufactured by:

Medi-Globe GmbH | Medi-Globe-Straße 1-5 | 83101 Achenmühle, Germany
☎ +49 80 32 / 973- 379 | 📠 +49 80 32 / 973-339 | ✉ sales@medi-globe.de