


## SCOPE TRANSPORT

Hygienic, efficient and easy to handle



## SCOPE:TRUCK

- » Shape of the tray is designed to dimensions for the endoscopes of all known manufacturers
- » Stable, plastic coated steel and aluminium construction
- » Four antistatic steering-wheels Ø 150 mm, two of them lockable
- » Trays made of stainless steel move on telescopic runners
- » SCOPE:TRUCK incl. 100 SCOPE:COVER (starter set)
- » (HxWxD): 1290 mm (with basin), 870 mm, 850 mm



7306-T00

## SCOPE:PORTER

The SCOPE:PORTER has been designed with modern Endoscopy Departments in mind:

- » Purpose-built with large basins for transporting six endoscopes
- » Reduces risk of contamination and scope damage by minimizing the handling of endoscopes during transport
- » Improved infection prevention because of better coverage from the large SCOPE:COVER (7306-COVER)
- » Easily maneuverable, even when fully loaded
- » Large, sturdy casters and quality steel and aluminum construction make narrow hallways and doors easy to navigate
- » Large work surface on the top of the SCOPE:PORTER makes it easy to load an endoscope
- » SCOPE:PORTER incl. 100 SCOPE:COVER (starter set)
- » (HxWxD): 1160 mm, 830 mm, 690 mm



7306-P00

## SCOPE:PORTER mini

- » Purpose-built with basins for transporting seven endoscopes
- » Reduces risk of contamination and scope damage by minimizing the handling of endoscopes during transport
- » Improved infection prevention because of better coverage from the small SCOPE:COVER (7316-COVER)
- » Easily maneuverable, even when fully loaded
- » Large, sturdy casters and quality steel and aluminum construction make narrow hallways and doors easy to navigate
- » Large work surface on the top of the SCOPE:PORTER makes it easy to load an endoscope
- » SCOPE:PORTER mini incl. 100 SCOPE:COVER (starter set)
- » (HxWxD): 1140 mm, 670 mm, 570 mm



7316-P00

# SCOPE:TRUCK & SCOPE:PORTER

7306-P00	SCOPE:PORTER for transport of six endoscopes incl. 100 SCOPE:COVER
7306-SET	2500 SCOPE:COVER (7306-COVER) incl. SCOPE:PORTER
7316-P00	SCOPE:PORTER mini for transport of seven endoscopes incl. 100 SCOPE:COVER
7316-SET	2500 SCOPE:COVER (7316-COVER) incl. SCOPE:PORTER mini
7306-T00	SCOPE:TRUCK for transport of six endoscopes incl. 100 SCOPE:COVER
7306-SET-Truck	2500 SCOPE:COVER (7306-COVER) incl. SCOPE:TRUCK



Trolleys



Endoscopy Accessories  
General Hospital Equipment



Electronic Accessories  
General Hospital Equipment



Flexible Accessories

ENDO-TECHNIK W. Griesat GmbH  
Heiligenstock 32  
42697 Solingen  
Germany

Tel.: +49 (0)212 / 267 43-0  
Fax: +49 (0)212 / 267 43-50

E-Mail: [info@endo-technik.de](mailto:info@endo-technik.de)  
Internet: [www.endo-technik.de](http://www.endo-technik.de)



2013.02

Whilst every effort is made to ensure that the specifications, descriptions and illustrations of all Products in this brochure are correct at the time of going to press, these are also subject to change without prior notice. Omissions excepted.