

Stone Extraction Baskets

Single Use

Stone Extraction Baskets

- 45° Injection Port & Integrated Sealing Cap
User-friendly and Easy Injection of Contrast Medium
- Detachable Instrument Holder & Outer Sheath
Applicable to Emergency Lithotripsy Spiral and StoneBuster Handle in Emergency Cases
- Ergonomic Instrument Holder
Simple and Safe Basket Manipulation



Stone Extraction Baskets

Single Use

Features & Benefits

- 45° Injection Port & Integrated Sealing Cap
User-friendly and Easy Injection of Contrast Medium
- Detachable Instrument Holder & Outer Sheath
Applicable to Emergency Lithotripsy Spiral and StoneBuster Handle in Emergency Cases
- Ergonomic Instrument Holder
Simple and Safe Basket Manipulation

Stone Extraction Basket Set – Set Consists of Basket, Outer Sheath (Working Length 180 cm), Sealing Cap and Instrument Holder

Catalog No.	Description	Basket Diameter	Outer Sheath Diameter
GEB-11-19-420	4-Wire Baskets	20 mm	1.9 mm
GEB-11-19-425	4-Wire Baskets	25 mm	1.9 mm
GEB-11-19-430	4-Wire Baskets	30 mm	1.9 mm
GEB-11-19-620	6-Wire Baskets	20 mm	1.9 mm
GEB-11-19-625	6-Wire Baskets	25 mm	1.9 mm
GEB-11-19-630	6-Wire Baskets	30 mm	1.9 mm
GEB-11-25-420	4-Wire Baskets	20 mm	2.5 mm
GEB-11-25-425	4-Wire Baskets	25 mm	2.5 mm
GEB-11-25-430	4-Wire Baskets	30 mm	2.5 mm
GEB-11-25-620	6-Wire Baskets	20 mm	2.5 mm
GEB-11-25-625	6-Wire Baskets	25 mm	2.5 mm
GEB-11-25-630	6-Wire Baskets	30 mm	2.5 mm

Ver. 4.0 (2010-08)

Distributed by:

Comercializadora Endosistemas Ltda.
Av. Suecia N° 3243, Ñuñoa, Stgo.
Fono: 22090262 - F/Fax: 27855826
ventas@endosistemas.cl

Manufactured By:

Medi-Globe GmbH Medi-Globe-Straße 1-5 83101 Achenmühle, Germany
☎ +49 80 32/973-379 📠 +49 80 32/973-339 ✉ sales@medi-globe.de
www.medi-globe.de