

Heyinovo

Endoscope Accessories

CE

Single use

Reusable



Made in China

Disposable Cleaning Brushes

Clean endoscope working channel and port hole.

Single use



Disposable Double Ended Cleaning Brushes



Channel cleaning brush.




Item No.	Brush Diameter & Working Length	Working Channel Diameter	Sterile
WS-1818WA	Ø 5.5 × 1800 (mm)	2.8~4.2 mm	
WS-1823WA	Ø 5.5 × 2300 (mm)	2.8~4.2 mm	
WS-1510WB	Ø 3.0 × 1000 (mm)	2.0~2.2 mm	
WS-1523WB	Ø 3.0 × 2300 (mm)	2.0~2.2 mm	
WS-0910WC	Ø 2.0 × 1000 (mm)	1.2 mm	
WS-0923WC	Ø 2.0 × 2300 (mm)	1.2 mm	




Item No.	Brush Diameter & Working Length	Working Channel Diameter	Sterile
WS-1818WG	Ø 5.5 × 1800 (mm)	2.8~4.2 mm	
WS-1823WG	Ø 5.5 × 2300 (mm)	2.8~4.2 mm	
WS-1510WH	Ø 3.0 × 1000 (mm)	2.0~2.2 mm	
WS-1523WH	Ø 3.0 × 2300 (mm)	2.0~2.2 mm	



Item No.	Brush Diameter & Working Length	Working Channel Diameter	Sterile
WS-1818WA(P)	Ø 5.5 × 1800 (mm)	2.8~4.2 mm	
WS-1823WA(P)	Ø 5.5 × 2300 (mm)	2.8~4.2 mm	
WS-1510WB(P)	Ø 3.0 × 1000 (mm)	2.0~2.2 mm	
WS-1518WB(P)	Ø 3.0 × 1800 (mm)	2.0~2.2 mm	
WS-1523WB(P)	Ø 3.0 × 2300 (mm)	2.0~2.2 mm	
WS-0910WC(P)	Ø 2.0 × 1000 (mm)	1.2 mm	
WS-0923WC(P)	Ø 2.0 × 2300 (mm)	1.2 mm	



Item No.	Brush Diameter & Working Length	Working Channel Diameter	Sterile
WS-1823WG/D	Ø 5.5/10 × 2300 (mm)	2.8~4.2 mm	



Item No.	Brush Diameter & Working Length	Working Channel Diameter	Sterile
WS-1823WA(P)/D	Ø 5.5/10 × 2300 (mm)	2.8~4.2 mm	



Item No.	Brush Diameter & Working Length	Working Channel Diameter	Sterile
WS-1523WA(P)/DH	Ø 5.5/10 × 2300 (mm)	2.0~4.2 mm	

Disposable Cleaning Brushes

Single use



Disposable Single Ended Cleaning Brushes



Channel cleaning brush.



Item No.	Brush Diameter & Working Length	Working Channel Diameter	Sterile
WS-1818VA	Ø 5.5 × 1800 (mm)	2.8~4.2 mm	
WS-1823VA	Ø 5.5 × 2300 (mm)	2.8~4.2 mm	
WS-1510VB	Ø 3.0 × 1000 (mm)	2.0~2.2 mm	
WS-1523VB	Ø 3.0 × 2300 (mm)	2.0~2.2 mm	
WS-0910VC	Ø 2.0 × 1000 (mm)	1.2 mm	
WS-0923VC	Ø 2.0 × 2300 (mm)	1.2 mm	



Item No.	Brush Diameter & Working Length	Working Channel Diameter	Sterile
WS-1818VG	Ø 5.5 × 1800 (mm)	2.8~4.2 mm	
WS-1823VG	Ø 5.5 × 2300 (mm)	2.8~4.2 mm	
WS-1510VH	Ø 3.0 × 1000 (mm)	2.0~2.2 mm	
WS-1523VH	Ø 3.0 × 2300 (mm)	2.0~2.2 mm	



Item No.	Brush Diameter & Working Length	Working Channel Diameter	Sterile
WS-1818VA(P)	Ø 5.5 × 1800 (mm)	2.8~4.2 mm	
WS-1823VA(P)	Ø 5.5 × 2300 (mm)	2.8~4.2 mm	
WS-1826VA(P)	Ø 5.5 × 2600 (mm)	2.8~4.2 mm	
WS-1510VB(P)	Ø 3.0 × 1000 (mm)	2.0~2.2 mm	
WS-1523VB(P)	Ø 3.0 × 2300 (mm)	2.0~2.2 mm	
WS-0910VC(P)	Ø 2.0 × 1000 (mm)	1.2 mm	
WS-0918VC(P)	Ø 2.0 × 1800 (mm)	1.2 mm	
WS-0923VC(P)	Ø 2.0 × 2300 (mm)	1.2 mm	
WS-0937VC(P)	Ø 2.0 × 3700 (mm)	1.2 mm	

Disposable Port Hole Cleaning Brushes



Port hole cleaning brush.



Item No.	Brush Diameter & Working Length	Sterile
WS-1001W1	Ø 10 × 60 (mm)	
WS-1001W2	Ø 10 × 60 (mm)	
WS-1001W8	Ø 5.5/10 × 40/60(mm)	

Distributed by:

ENDO SISTEMAS
El valor de la experiencia

Comercializadora Endosistemas Ltda.

Av. Suecia Nº 3243, Ñuñoa, Santiago

Tel.: 22090262 - Fax.: 27855826

ventas@endosistemas.cl

www.endosistemas.cl