

Heyinovo

Endoscope Accessories

CE

Single use

Reusable



Disposable Biopsy Forceps

Collect tissue from human body through endoscope.



Disposable Oval Cups Biopsy Forceps



Oval, fenestrated cups for deeper biopsy.

Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Sterile
WS-1810BT	Ø 1.8 × 1000 (mm)	Ø 2.0 (mm)	●
WS-1815BT	Ø 1.8 × 1500 (mm)	Ø 2.0 (mm)	●
WS-2415BT	Ø 2.4 × 1500 (mm)	Ø 2.8 (mm)	●
WS-2422BT	Ø 2.4 × 2200 (mm)	Ø 2.8 (mm)	●

Disposable Oval Cups With Needle Biopsy Forceps



Oval, fenestrated cups with needle for accurate targeting, no sliding, for gathering bigger tissue sample.

Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Sterile
WS-2415BW	Ø 2.4 × 1500 (mm)	Ø 2.8 (mm)	●
WS-2422BW	Ø 2.4 × 2200 (mm)	Ø 2.8 (mm)	●

Disposable Alligator Cups Biopsy Forceps



Serrated, Alligator cups for deep biopsy of hard and difficult to access tissue.

Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Sterile
WS-2415BU	Ø 2.4 × 1500 (mm)	Ø 2.8 (mm)	●
WS-2422BU	Ø 2.4 × 2200 (mm)	Ø 2.8 (mm)	●

Disposable Alligator Cups With Needle Biopsy Forceps



Serrated, Alligator cups for deep biopsy of hard and difficult to access tissue.

With needle for accurate targeting, no slippage.

Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Sterile
WS-2415BN	Ø 2.4 × 1500 (mm)	Ø 2.8 (mm)	●
WS-2422BN	Ø 2.4 × 2200 (mm)	Ø 2.8 (mm)	●

Biopsy Forceps (Reusable)

Collect tissue from human body through endoscope.



Standard Biopsy Forceps

Autoclavable 134°C

Standard, round and fenestrated cups for minimal tissue damage.



Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Autoclavable
WF-1807BS	Ø 1.8× 700 (mm)	Ø 2.0 (mm)	●
WF-1810BS	Ø 1.8× 1000 (mm)	Ø 2.0 (mm)	●
WF-1815BS	Ø 1.8× 1500 (mm)	Ø 2.0 (mm)	●
WF-1825BS	Ø 1.8× 2500 (mm)	Ø 2.0 (mm)	●
WF-2415BS	Ø 2.4× 1500 (mm)	Ø 2.8 (mm)	●
WF-2423BS	Ø 2.4× 2300 (mm)	Ø 2.8 (mm)	●

* Other specification can be customized according to requirement.

Oval Fenestrated Cups

Autoclavable 134°C

Oval, fenestrated cups for deeper biopsy.



Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Autoclavable
WF-1807BT	Ø 1.8× 700 (mm)	Ø 2.0 (mm)	●
WF-1810BT	Ø 1.8× 1000 (mm)	Ø 2.0 (mm)	●
WF-1815BT	Ø 1.8× 1500 (mm)	Ø 2.0 (mm)	●
WF-1825BT	Ø 1.8× 2500 (mm)	Ø 2.0 (mm)	●
WF-2415BT	Ø 2.4× 1500 (mm)	Ø 2.8 (mm)	●
WF-2423BT	Ø 2.4× 2300 (mm)	Ø 2.8 (mm)	●

* Other specification can be customized according to requirement.

Oval Fenestrated Cups With Needle

Autoclavable 134°C

Oval, fenestrated cups with needle for accurate targeting, no sliding. Desinged to minimize tissue damage.



Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Autoclavable
WF-2415BW	Ø 2.4× 1500 (mm)	Ø 2.8 (mm)	●
WF-2423BW	Ø 2.4× 2300 (mm)	Ø 2.8 (mm)	●

* Other specification can be customized according to requirement.

Biopsy Forceps (Reusable)

Collect tissue from human body through endoscope.



Alligator Jaws Type Biopsy Forceps

Autoclavable 134°C

Standard, alligator jaws cups for tissue collection on timour.



Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Autoclavable
WF-1807BU	Ø 1.8 × 700 (mm)	Ø 2.0 (mm)	●
WF-1810BU	Ø 1.8 × 1000 (mm)	Ø 2.0 (mm)	●
WF-1815BU	Ø 1.8 × 1500 (mm)	Ø 2.0 (mm)	●
WF-1825BU	Ø 1.8 × 2500 (mm)	Ø 2.0 (mm)	●
WF-2415BU	Ø 2.4 × 1500 (mm)	Ø 2.8 (mm)	●
WF-2423BU	Ø 2.4 × 2300 (mm)	Ø 2.8 (mm)	●

* Other specification can be customized according to requirement.

Alligator, Single Moving Cup Biopsy Forceps

Autoclavable 134°C

Single moving cup for wall biopsies and biopsies of depressed lesions.



Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Autoclavable
WF-1807BV	Ø 1.8 × 700 (mm)	Ø 2.0 (mm)	●
WF-1810BV	Ø 1.8 × 1000 (mm)	Ø 2.0 (mm)	●
WF-1815BV	Ø 1.8 × 1500 (mm)	Ø 2.0 (mm)	●
WF-1825BV	Ø 1.8 × 2500 (mm)	Ø 2.0 (mm)	●
WF-2415BV	Ø 2.4 × 1500 (mm)	Ø 2.8 (mm)	●
WF-2423BV	Ø 2.4 × 2300 (mm)	Ø 2.8 (mm)	●

* Other specification can be customized according to requirement.

Manufactured by:



Wilson Instrument (SHA) Co., Ltd.
25D, Heyi Bussines Plaza N° 420
Jiang Ning Rd, Shangai, RP China
Tel: ++86 21 52930848
Fax: ++86 21 52930146
heyinovo@gmail.com
www.heyinovo.com

Distributed by:

Comercializadora Endosistemas Ltda.
Suecia N° 3243, Ñuñoa, Santiago
Tel.: 22090262 - Fax.: 27855826
ventas@endosistemas.cl
www.endosistemas.cl