

Heyinovo

Endoscope Accessories

CE

Single use

Reusable



Foreign Body Removal (Reusable)

Grasp foreign body through endoscope.



Rat Tooth Type Grasping Forceps

Autoclavable 134°C

Rat tooth for flat and hard objects (e.g. Coins).



Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Autoclavable
WF-1807GI	Ø 1.8 × 700 (mm)	Ø 2.0 (mm)	●
WF-1810GI	Ø 1.8 × 1000 (mm)	Ø 2.0 (mm)	●
WF-1815GI	Ø 1.8 × 1500 (mm)	Ø 2.0 (mm)	●
WF-1825GI	Ø 1.8 × 2500 (mm)	Ø 2.0 (mm)	●
WF-2415GI	Ø 2.4 × 1500 (mm)	Ø 2.8 (mm)	●
WF-2423GI	Ø 2.4 × 2300 (mm)	Ø 2.8 (mm)	●

Rat Tooth With Alligator Jaws Type Grasping Forceps

Autoclavable 134°C

For flat and hard objects, sharp tooth, long jaws for difficult retrieval.



Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Autoclavable
WF-1807GJ	Ø 1.8 × 700 (mm)	Ø 2.0 (mm)	●
WF-1810GJ	Ø 1.8 × 1000 (mm)	Ø 2.0 (mm)	●
WF-1815GJ	Ø 1.8 × 1500 (mm)	Ø 2.0 (mm)	●
WF-1825GJ	Ø 1.8 × 2500 (mm)	Ø 2.0 (mm)	●
WF-2415GJ	Ø 2.4 × 1500 (mm)	Ø 2.8 (mm)	●
WF-2423GJ	Ø 2.4 × 2300 (mm)	Ø 2.8 (mm)	●

Alligator Type Grasping Forceps

Autoclavable 134°C

Alligator jaws for sharp edged objects.



Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Autoclavable
WF-1807GH	Ø 1.8 × 700 (mm)	Ø 2.0 (mm)	●
WF-1810GH	Ø 1.8 × 1000 (mm)	Ø 2.0 (mm)	●
WF-1815GH	Ø 1.8 × 1500 (mm)	Ø 2.0 (mm)	●
WF-1825GH	Ø 1.8 × 2500 (mm)	Ø 2.0 (mm)	●
WF-2415GH	Ø 2.4 × 1500 (mm)	Ø 2.8 (mm)	●
WF-2423GH	Ø 2.4 × 2300 (mm)	Ø 2.8 (mm)	●

Foreign Body Removal (Reusable)

Grasp foreign body through endoscope.



Pelican Type Grasping Forceps

Autoclavable **134°C**

For big object. Opening width: 15 mm.



Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Autoclavable
WF-1807GK	Ø 1.8× 700 (mm)	Ø 2.0 (mm)	●
WF-1810GK	Ø 1.8× 1000 (mm)	Ø 2.0 (mm)	●
WF-1815GK	Ø 1.8× 1500 (mm)	Ø 2.0 (mm)	●
WF-1825GK	Ø 1.8× 2500 (mm)	Ø 2.0 (mm)	●
WF-2415GK	Ø 2.4× 1500 (mm)	Ø 2.8 (mm)	●
WF-2423GK	Ø 2.4× 2300 (mm)	Ø 2.8 (mm)	●

Rubber Tips Type Grasping Forceps

Autoclavable **134°C**

Rubber tips for grasping slippery, flat or slightness objects. (E.g. needles slightness wires).



Item No.	Max. Insertion Portion Diameter & Working Length	Working Channel Diameter	Autoclavable
WF-1807GN	Ø 1.8× 700 (mm)	Ø 2.0 (mm)	●
WF-1810GN	Ø 1.8× 1000 (mm)	Ø 2.0 (mm)	●
WF-1815GN	Ø 1.8× 1500 (mm)	Ø 2.0 (mm)	●
WF-1825GN	Ø 1.8× 2500 (mm)	Ø 2.0 (mm)	●
WF-2415GN	Ø 2.4× 1500 (mm)	Ø 2.8 (mm)	●
WF-2423GN	Ø 2.4× 2300 (mm)	Ø 2.8 (mm)	●

Manufactured by:



Wilson Instrument (SHA) Co., Ltd.
25D, Heyi Bussines Plaza N° 420
Jiang Ning Rd, Shangai, RP China
Tel: ++86 21 52930848
Fax: ++86 21 52930146
heyinovo@gmail.com
www.heyinovo.com

Distributed by:



Comercializadora Endosistemas Ltda.
Av. Suecia N° 3243, Ñuñoa, Santiago
Tel.: 22090262 - Fax.: 27855826
ventas@endosistemas.cl
www.endosistemas.cl